ENGLISH

**RANGE OF MOTION ADJUSTMENTS:**

The Catalyst•Elite comes with standard 5 degree extension stops. Additional flexion and extension stops may be installed if they are prescribed by your physician or needed for optimal performance or fit of the brace. These should be inserted by the technician who fits your brace. Optional 0°, 10°, 20°, 30°, and 40° extension stops and optional 45°, 60°, 75° and 90° flexion stops are included with each brace.

**IMPORTANT: Both hinges must be set to same extension and/or flexion settings. For 0° extension limiting, 0° extension stops must be used.**

**TO CHANGE AN EXTENSION STOP:**

Unscrew the screw on the front of the hinge cap. Shake out the 5 degree extension stop. Insert another extension stop into the front of the hinge and replace the screw. Repeat for the opposite hinge.

**TO INSERT A FLEXION STOP:**

Insert a flexion stop into the back of the hinge and use the extra screw that is attached to the flexion/extension stop tree to secure. Repeat for the opposite hinge.

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

ENGLISH

**INDICATIONS:**

- Stabilization of ACL, MCL, PCL, LCL or treatment of combined ligament instabilities
- Protection for surgical repairs
- Conservative treatment of ligament damage or deficiency
- Prophylactic bracing for athletic activities

**WARNINGS:**

If you experience pain, swelling, sensation changes, or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at 800.452.7993 or consult a physician.

**FOR SINGLE PATIENT USE ONLY.**

**NOTICE:** While every effort has been made using state-of-the-art techniques to obtain maximum compatibility of function, strength, durability and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through use of this device.

**CAUTION:** Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

**WARRANTY:** VQ OrthoCare warranties all its products from the original date of purchase against defects in materials and workmanship. Normal wear and tear during use of a product is not considered a defect. Contact VQ OrthoCare at **949.261.3000** for specific product warranty information.

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MANUFACTURED AT  
1390 Decision St, Suite A  
Vista, CA 92081, USA

© 2015 VisionQuest Industries Inc.  
VQ0151706REVD

# CATALYST•ELITE APPLICATION INSTRUCTIONS

(Instrucciones De Aplicación De Catalyst•Elite / Instructions d'application Catalyst•Elite)

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## ENGLISH

### STEP 1. UNFASTEN STRAPS

- To simplify application, attach each strap back onto itself to minimize tangling (Fig. 1).

### STEP 2. POSITION BRACE

- While sitting at edge of a chair, step through the brace and flex the knee to approximately 45°.

- Align hinge centers with top of kneecap and position slightly behind midline of leg (Fig. 2).

- Catalyst•Elite braces with the Active Thigh Cuff (ATC) may be adjusted further if the uprights are not positioned along the midline of the thigh or if there is gapping between the soft tissue and the ATC straps.

**NOTE:** ATC adjustment should be done by your healthcare professional or brace technician. Remove the screws for each strap. Adjust as needed to fit flush against the thigh with the side flanges sitting mid-medial and mid-lateral or slightly posterior. Replace screws. Trim as needed.

### STEP 3. FASTEN STRAPS

- While maintaining the brace in line with the top of the patella, fasten straps 1-5 in sequence indicated on strap tabs and in Fig. 3. Straps should be comfortably tightened to properly secure brace. Do not overtighten.

① Fasten lowest shin strap securely.

② Fasten upper thigh strap securely.

③ Fasten lower thigh strap securely.

④ Fasten upper shin strap securely.

⑤ Fasten calf strap securely. This strap is most important in suspending the brace. Make sure that it sits on top of the calf muscle.

**NOTE:** For PCL configuration (Fig. 4a), steps 1 & 2 are the same. The strap on front of thigh will be fastened 3rd. The upper shin strap is fastened 4th and the posterior calf strap is fastened last.

**NOTE:** For Combined Instability configuration (Fig. 4b), steps 1-5 are the same. The strap on the front of the thigh will be fastened last.

**NOTE:** D-ring liner tabs may be trimmed for comfort.

### STEP 4. CHECK FIT

- Stand up and walk around in the brace.

- Check to make sure that hinge centers are aligned with the middle of kneecap or slightly above and positioned slightly behind midline of leg. It is normal for the brace to settle into position. Make adjustments as necessary to correct brace position and ensure a comfortable fit.

### ATC MODELS

If desired, install ATC pivot clip where the ATC components cross, as shown. The pivot clip allows the ATC to pivot in the center during use and is there as an option if the crossing straps are moving too much or shift out of their proper place on the leg.

Also included is the rigid ATC straps that can be used in place of the flexible straps already installed if additional rigidity is desired.

### CARE AND MAINTENANCE

- Salt water, sand, dirt and other elements can cause damage. Thoroughly rinse brace with fresh water after exposure to these conditions. Towel or air dry.

- Liners should be removed and hand washed as needed using cold water and mild detergent. Rinse thoroughly and air dry (Do not machine wash. Do not machine dry).

- Lubricate hinge with dry silicone spray as needed

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