## ENGLISH

### INDICATIONS:

- Mild to moderate unicompartmental knee osteoarthritis with ligament instabilities
- Increased rehabilitative and functional support for diagnoses that require reduced medial or lateral joint loads

### WARNINGS:

If you experience pain, swelling, sensation changes, or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at **800.452.7993** or contact the provider that delivered your brace or consult a physician.

### FOR SINGLE PATIENT USE ONLY.

**NOTICE:** While every effort has been made using state-of-the-art techniques to obtain maximum compatibility of function, strength, durability and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through use of this device.

**CAUTION:** Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

**WARRANTY:** VQ OrthoCare warranties all its products from the original date of purchase against defects in materials and workmanship. Normal wear and tear during use of a product is not considered a defect. Contact VQ OrthoCare at **949.261.3000** for specific product warranty information. Outside the US, please contact the provider that delivered your brace.

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**MISE EN GARDE:** Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

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Subject to one or more patents listed on  
www.vqorthocare.com/patents

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MANUFACTURED AT  
1390 Decision St, Suite A  
Vista, CA 92081, USA

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



## FreeFlex OA

# FREEFLEX OA APPLICATION INSTRUCTIONS

(Instrucciones de Aplicación de FreeFlex OA / Instructions de Mise en Place de FreeFlex)

## ENGLISH

### STEP 1. POSITION BRACE

- Sit with knee **slightly** (5 - 10°) bent.
- Position the brace by placing the hinge on the side of the knee **opposite** of your affected compartment.

**Medial OA:** Place hinge on **OUTSIDE** of knee  
**Lateral OA:** Place hinge on **INSIDE** of knee

**NOTE:** The brace needs to be applied directly to skin and should **NOT** be worn over clothing.

- Hold the brace by the hinge and position the hinge center (marked by the grey extension stop) even with the **TOP** of your kneecap and slightly posterior (towards the back) to your leg.

### STEP 2. FASTEN WRAP

- Secure the top and bottom wrap on your leg according to the numbers on the inside of the brace, maintaining the hinge in position. (See **Fig. 1 & Fig. 2**)

- Secure the middle two straps according to the numbers on the inside of the brace above and below your knee cap. (See **Fig. 1 & Fig. 3**)

### STEP 3. FASTEN STRAPS

- Pull the top strap behind the knee and secure to the **BOTTOM** buckle. Tighten strap using Velcro pull. (See **Fig. 4**) Do not overtighten.

- Pull the bottom strap behind the knee and secure to the **TOP** buckle. Tighten strap using Velcro pull. (See **Fig. 5**) Do not overtighten. Loosen straps if you feel pinching.

- Brace should look like **Fig. 6** from behind the knee.

### APPLY CORRECTION

The criss-cross straps provide some OA correction. The Q-hinge will provide additional correction. You should only feel gentle pressure – too much correction can cause pain. Insert the allen wrench

into the Q-hinge and turn either clockwise or counter-clockwise until the hinge has the same shape or curve as your leg (neutral). There should be no pressure at all on the knee at this time. Walk around and make sure the brace is comfortable and that you don't feel pressure.

Next, sit back down and bend knee at least 30-degrees, then insert the allen wrench again into the Q-hinge and turn away from your body about a 1/4 turn (90°). This is approximately 3 degrees of correction and you should feel gentle pressure during full-extension. If you do not feel **ANY** pressure, try another 1/8 to 1/4 turn. Stand up and walk around in the brace. Do **NOT** overcorrect as this will cause discomfort and may increase your pain over time.

### “WEAR TIME”

Gradually increase brace wear time to allow your body to adjust to the brace. Start with 30 minutes on day one and gradually increase. Pain caused by osteoarthritis may take time to subside. It is important that you routinely wear the brace when you're on your feet as it will reduce joint loading and allow you to begin feeling relief. It may take a number of weeks for relief to become significant.

### CARE AND CLEANING

Gently hand wash material in cold water with mild detergent. Air dry. Do **NOT** machine wash or dry.

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© 2015 VisionQuest Industries, Inc.