Home Exercise Kit

Foot and Ankle

User Manual

VQ OrthoCare Exercise Kit

Foot and Ankle Kit


Items included in your kit:

- Marbles

- Pro-Stretch®

- Jump Rope

- **Exercise Band Kit**
 - 2 Handles
 - Blue, Green, Red Rubber Bands
 - Strap with Loop
 - Pouch

- **Towel**

VQ OrthoCare is not liable for any misuse or misunderstanding of the exercise product or operating manual. Please call your local representative if any additional assistance is required regarding this product and its operating instructions.

Caution: Consult with your therapist or physician before starting any rehabilitation or exercise program. The products in this exercise kit are not to be used by children without adult supervision.

Caution: This product contains natural rubber latex which may cause allergic reactions. As with all natural rubber products, keep this product away from grease, oil or excessive heat.

Disclaimer: Equipment and supplies illustrated and specified are subject to modification consistent with product availability or other factors.

VQ OrthoCare

18011 Mitchell South, Irvine, CA 92614-6007

800.266.6969 • Fax 800.821.8012

www.vqorthocare.com • customerservice@vqorthocare.com

The publishers acknowledge that the illustrations used in this publication on pages 2,3,4,6,11,12,13,14,15,16,17 are the copyrighted property of VHI Inc.


Pro-Stretch is a registered trademark of Prism Enterprises.

This operating manual is proprietary and confidential. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording or information storage and retrieval system, without permission in writing from VQ OrthoCare.

Quick Reference Exercise Listing


Foot and Ankle

Ankle Pumps	2
Ankle Alphabet	3
Towel Gather	4
Marble Pick Up	5
Foot Stretches	6
Stretch Exercise – Calf	7
Stretch Exercise – Soleus	8
Plantarflexion/Dorsiflexion	9
Inversion/Eversion	10
Heel Raises	11
Toe Raises	12
Single Leg Balance	13
Heel/Toe Walking	14
Exercise Band Exercises – Dorsiflexion	15
Exercise Band Exercises – Plantarflexion	16
Exercise Band Exercises – Inversion	17
Jump Rope	18


ANKLE PUMPS

Sit or lie with legs stretched out in front of you. Toes should be pointing up. One foot at a time, first point your toes down and away from you, then flex your ankle bring toes toward you as far as comfortable. Repeat exercise for a total of 10 repetitions. Complete 3 sets of 10 repetitions each. Repeat entire exercise with the opposite foot. Do this exercise 3-4 times per day.


ANKLE ALPHABET

Sit with legs freely hanging off the end of a table or lay on the floor with legs stretched out in front of you. Point your toes down and away from you. Move toes and ankle as if you were writing the alphabet. This will move your ankle in all directions. Repeat exercise twice. Do this exercise 3-4 times per day.


TOWEL GATHER

Sit in a chair with your feet flat on the floor. Spread a towel on the floor. Gather the towel with your toes until the towel is fully gathered around your foot. Repeat 10 times. Do this exercise 3 times per day.


MARBLE PICK UP

Sit in a chair with your feet flat on the floor. Place marbles in a pile in front of you. Pick up a few marbles between your toes and create a separate pile of marbles a few inches away from the first pile. Continue until all the marbles have been transferred to this new pile. Rebuild the pile 5 times. Do this exercise 3 times daily.


FOOT STRETCHES

Sit in a chair or on a firm surface with your injured leg stretched out in front of you. Loop a towel or exercise band around your foot and pull towards you until you feel a gentle stretch. Hold the stretch for 20 seconds and relax for 20 seconds. Next pull the towel or exercise band at an angle so your foot is turned to the inside. Hold the stretch for 20 seconds and relax for 20 seconds. Repeat the exercise pulling the towel or exercise band in the opposite direction so that the foot is turned to the outside. Finally, bend your injured leg to 90°, loop the towel or exercise band behind your toes and pull towards you. Hold the stretch for 20 seconds and relax for 20 seconds. Do this exercise 3 times per day.


STRETCH EXERCISE – CALF

Place the stretching device perpendicular to and an arm's length from a wall or other supportive surface. Face the wall and place both hands against the wall. Step onto the stretching device straight up with all of your weight and lift the opposite foot off the floor keeping your shoulders, hips and ankle in line. Keeping the leg in the stretching device straight, rock back and forth until you feel a slight tension or stretch in your upper calf muscle. Hold for 30 seconds and relax for 5 seconds. Complete 3 sets of 30 second stretches each. Switch to the other leg and repeat. Do this exercise 3 times daily.


STRETCH EXERCISE – SOLEUS

Position your foot in the stretching device following the instructions under "Stretch Exercise – Calf." Once you feel a slight stretch in your calf muscle, bend your knee slightly keeping your heel back. Hold for 30 seconds and relax for 5 seconds. Complete 3 sets of 30 second stretches each. Switch to the other leg and repeat. Do this exercise 3 times daily.


PLANTARFLEXION / DORSIFLEXION

Sit in a chair placing your foot into the stretching device. Slowly move your ankle forward and push your toes and forefoot toward the floor. Then rock back pushing your heel toward the floor. Repeat this motion 20 times. Complete 3 sets of 20 repetitions each. Do this exercise 3 times a day. NOTE: You can progress to performing this exercise while standing with most of your weight on your uninjured leg. Slowly put more weight on the injured leg as you become stronger.

INVERSION / EVERSION

Sit in a chair, place the stretching device sideways in front of you. Place the middle of your foot on the stretching device and rock your foot side to side. Repeat this motion 20 times. Complete 3 sets of 20 repetitions each. Do this exercise 3 times a day. NOTE: You can progress to performing this exercise while standing, always using a stable surface for support.


HEEL RAISES

Stand behind a chair and balance yourself with your feet slightly apart. Rise up on your toes, hold for 5 seconds and relax. Repeat 10 times. Complete 3 sets of 10 repetitions each. Repeat this exercise with your toes pointed inward then with your toes pointed outward. Do this exercise 3 times daily.

TOE RAISES

Stand with your back approximately 6 inches from a wall. Rock back on your heels so that your toes rise slightly off the ground without leaning back and touching the wall for support. Hold for 5 seconds. Repeat 10 times. Complete 3 sets of 10 repetitions each. Do this exercise 3 times daily.


SINGLE LEG BALANCE

Stand behind a chair and balance yourself on your uninjured leg, trying not to use the chair for balance or support. Begin the exercise with your eyes open, progressing to eyes closed and bent knee position. Hold the single leg balance position for 30 seconds. Repeat 3 times. Do this exercise 3-4 times daily.


HEEL/TOE WALKING

Walk on heels only for a total of 30 steps. Walk on toes only for a total of 30 steps. Do this exercise 3 times daily.

EXERCISE BAND – DORSIFLEXION

Anchor the exercise band into the door by slipping the band between the door and the frame approximately 6 inches from the floor and closing the door. Sit on the floor with your legs straight out in front of you, toes pointed towards the door. Loop the strap around the top of your foot. Keeping your knee straight, pull your toes toward you. Slowly return to the starting position. Repeat 10 times. Complete 3 sets of 10 repetitions each. Do this exercise 3 times daily.

EXERCISE BAND – PLANTARFLEXION

Sit on the floor with your legs straight out in front of you. Loop one end of the exercise band around the top of your foot just below the toes. Hold the other end of the band firmly in both hands. Gently press the top of your foot down, pointing your toes away from you and stretching the band. Return to the starting position. Complete 3 sets of 10 repetitions each. Do this exercise 3 times daily.

EXERCISE BAND INVERSION

Sit on the floor with your legs straight out in front of you. Cross your uninjured leg over your injured ankle. Slip the top of your foot on the injured side through the exercise band and then loop the band around your uninjured foot so the band is anchored. Hold the other end of the band firmly in your hands. Turn your injured foot inward and upward, stretching the band. Return to the starting position. Complete 3 sets of 10 repetitions each. Do this exercise 3 times daily.

JUMP ROPE

Jump rope landing on both legs. Try to build up your endurance to 5 minutes. Note: As you build strength, progress to landing only on the injured leg jumping for up to 5 minutes.

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