Home Exercise Kit

Knee

User Manual

VQ OrthoCare Exercise Kit

Knee Kit

Items included in your kit:

- **Air Pillow**

- **Water Weight Bag**

- **Exercise Band Kit**
 - 2 Handles
 - Blue, Green, Red Rubber Bands
 - Strap with Loop
 - Pouch

- **Towel**

VQ OrthoCare is not liable for any misuse or misunderstanding of the exercise product or operating manual. Please call your local representative if any additional assistance is required regarding this product and its operating instructions.

Caution: Consult with your therapist or physician before starting any rehabilitation or exercise program. The products in this exercise kit are not to be used by children without adult supervision.

Caution: This product contains natural rubber latex which may cause allergic reactions. As with all natural rubber products, keep this product away from grease, oil or excessive heat.

Disclaimer: Equipment and supplies illustrated and specified are subject to modification consistent with product availability or other factors.

VQ OrthoCare

18011 Mitchell South, Irvine, CA 92614-6007

800.266.6969 • Fax 800.821.8012

www.vqorthocare.com • customerservice@vqorthocare.com


The publishers acknowledge that the illustrations used in this publication on pages 2,3,4,5,6,7,10, 11,12,13,15,16,17 are the copyrighted property of VHI Inc.

This operating manual is proprietary and confidential. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording or information storage and retrieval system, without permission in writing from VQ OrthoCare.

Quick Reference Exercise Listing


Knee

Quad Sets	2
Ankle Pumps	3
Hamstring Sets	4
Elevated Knee Extensions	5
Gluteal Sets	6
Knee Flexion	7
Heel Slides	8
Hamstring Stretches	9
Heel/Calf Stretches with Band	10
Weighted Exercises	
Straight Leg Raise	11
Knee Extensions	12
Hamstring Curls	13
Abduction	14
Adduction	15
Exercise Band Exercises	
Seated Knee Extensions	16
Seated Hamstring Curls	17


QUAD SETS

Sit on a flat surface with your legs straight out in front of you. Tighten the thigh muscle of your injured leg by pushing the back of your knee down into the flat surface. Hold for a count of 5 seconds. Relax. Do 20 repetitions. You may use the folded towel under your ankle for support when performing this exercise. Repeat this exercise 3 times a day.


ANKLE PUMPS

Sit on a flat surface with your leg straight out in front of you. Bend your ankle up and down by pointing your toes and then pulling your toes back toward you. Do 20 repetitions and repeat 3 times a day.


HAMSTRING SETS

Sit on a flat surface with your uninjured leg straight out in front of you and your injured leg bent slightly at the knee. Push the heel of your injured foot into the floor. You may use the folded towel to support your heel. You may also place the inflated roll under your injured knee to support it while doing this exercise. Hold for a count of 5 seconds. Relax. Do 20 repetitions 3 times a day.

ELEVATED KNEE EXTENSIONS

Lie on your back with your legs straight out in front of you. Place the air pillow under your injured knee. Straighten your knee by tightening your thigh muscle and hold for a count of 5 seconds. Relax. Do not lift your leg off the roll. Do 20 repetitions 3 times a day.


GLUTEAL SETS

Stand with your feet together. Squeeze your buttocks together and hold for a count of 5 seconds. Relax. Do 20 repetitions 3 times a day.


KNEE FLEXION

Sit in a chair and lean back against the backrest. Use the blue exercise band without the handles attached. Place the band over the top of your foot at the ankle. Hold the ends of the band with both hands and pull on the band slowly to bend your injured knee and slide your foot as far back under the chair as possible. Hold for 3 seconds. Relax. **DO NOT LEAN FORWARD** while performing this exercise. Do 20 repetitions 3 times a day.


HEEL SLIDES

Lie flat on your back with your injured knee bent and foot flat on the floor. Slowly slide your heel toward your buttock while keeping your knee pointed up at the ceiling. Slowly slide your leg back down to a straight position. Do 20 repetitions 3 times a day.


HAMSTRING STRETCHES

Sit with one leg off the bed with the foot flat on the floor. Keep the other leg straight on the bed. Lean forward reaching for your toes of the straightened leg. Keep your back straight. Hold for a count of 15-20 then relax. Do not bounce forward. Alternate and stretch opposite leg. Do 10 repetitions 3 times a day.

HEEL/CALF STRETCHES WITH BAND

Sit on a bed or floor with your legs out in front of you. Start with the uninjured leg then progress to the injured leg. Use the blue exercise band without the handles attached. Place the band around the ball of the forefoot. Hold the ends of the band with both hands and pull on the band slowly until you feel a good stretch (you should feel a stretch in your heel, calf and behind the knee). Hold the stretch for 15 - 20 seconds then relax. Do 10 repetitions 3 times a day.

STRAIGHT LEG RAISE

Lie flat on your back with your legs extended and toes pulled back. Keeping your injured leg straight, slowly lift your leg up by tightening your thigh muscle. Make sure to lift the heel first. Lift leg about 6 inches off the floor. Hold the leg in the raised position for 3 seconds. Slowly lower your leg back down. Relax. Do 10-20 repetitions 3 times a day.

Optional Water Weight Straight Leg Raise:

Fill both sides of the weight bag with equal amounts of water. Start at 1 pound on each side and increase amount as you increase strength. Put your forefoot through the hole in the weight bag and drape the water bag around the foot of your injured leg. Do exercise as instructed above.


KNEE EXTENSIONS

Sit in a chair or at the end of a table so that your leg is hanging above the floor. Cross the ankle of your uninjured leg over the ankle of the injured leg. Straighten your injured leg, slowly lift your leg up by tightening your thigh muscle. Hold for 3 seconds then lower your leg slowly. Relax. Do 3 sets, 10 repetitions.

Optional Water Weight Knee Extensions:

Fill both sides of the weight bag with equal amounts of water. Start at 1 pound on each side and increase amount as you increase strength. Put your forefoot through the hole in the weight bag and drape the water bag around the foot of your injured leg. Do exercise as instructed above.


HAMSTRING CURLS

Lie face down. Bend your injured knee. Squeeze your muscle in the back part of your leg bringing your heel toward your buttocks and hold for 3 seconds. Slowly lower your leg. Relax. Do 3 sets, 10 repetitions.

Optional Water Weight Hamstring Curls:

Fill both sides of the weight bag with equal amounts of water. Start at 1 pound on each side and increase amount as you increase strength. In sitting position place your forefoot through the hole in the weight. Turn over and drape the water weight bag over the sole and heel of your foot. Do exercise as instructed above.


ABDUCTION

Lay on your side with your injured leg on top. Keep your uninjured leg straight and raise your injured leg up slowly as far as you can go. Hold for 3 seconds. Slowly lower your leg. Relax. Do 3 sets, 10 repetitions.

Optional Water Weight Abduction Exercises:

Fill both sides of the weight bag with equal amounts of water. Start at 1 pound on each side and increase amount as you increase strength. Put your forefoot through the hole in the weight bag and drape the water bag around the foot and ankle of your injured leg. Do exercise as instructed above.


ADDUCTION

Lay on your side with your uninjured leg on top. Cross your uninjured leg over the top of your injured leg and place your foot flat on the floor. Raise your uninjured leg as far as you can go while keeping the leg straight. Hold for 3 seconds. Slowly lower your leg. Relax. Do 10 repetitions 3 times a day.


Optional Water Weight Abduction Exercises:

Fill both sides of the weight bag with equal amounts of water. Start at 1 pound on each side and increase amount as you increase strength. Put your forefoot through the hole in the weight bag and drape the water bag **around the foot and** ankle of your injured leg. Do exercise as instructed above.

SEATED KNEE EXTENSIONS

Replace the exercise band handle with the exercise band strap provided. Place the handle end between the door and the doorframe about 2 inches from the floor and close the door to secure the band. Place the exercise band strap around the ankle of the injured leg. Sit in a chair with your back to the door far enough away so the strap is suspended under the chair, above the ground and slightly taut. Place your hands at your sides or grip the chair. Place your uninjured foot flat on the ground. Squeeze your thigh and slowly straighten the leg with the band so that your toes point toward the ceiling. Hold for 5 seconds. Slowly lower your foot back to the floor. Relax. Do 3 sets, 10 repetitions.

SEATED HAMSTRING CURLS

Replace the exercise band handle with the exercise band strap provided. Place the handle end between the door and the doorframe about 2 inches from the floor and close the door to secure the band. Place the exercise band strap around the ankle of the injured leg. Sit in a chair facing the door with your uninjured foot flat on the ground. With your injured leg in the exercise band strap and extended toward the door the band should be taut. Slowly pull your leg back underneath the chair as far as possible. Hold for 5 seconds and slowly return leg to the original position. Relax. Do 3 sets, 10 repetitions.

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