Home Exercise Kit

Shoulder

User Manual

VQ OrthoCare Exercise Kit

Shoulder Kit

Items included in your kit:

- Wand

- Pulley with Handles

- **Exercise Band Kit**
 - 2 Handles
 - Blue, Green, Red Rubber Bands
 - Strap with Loop
 - Pouch

- Towel

VQ OrthoCare is not liable for any misuse or misunderstanding of the exercise product or operating manual. Please call your local representative if any additional assistance is required regarding this product and its operating instructions.

Caution: Consult with your therapist or physician before starting any rehabilitation or exercise program. The products in this exercise kit are not to be used by children without adult supervision.

Caution: This product contains natural rubber latex which may cause allergic reactions. As with all natural rubber products, keep this product away from grease, oil or excessive heat.

Disclaimer: Equipment and supplies illustrated and specified are subject to modification consistent with product availability or other factors.

VQ OrthoCare

18011 Mitchell South, Irvine, CA 92614-6007

800.266.6969 • Fax 800.821.8012

www.vqorthocare.com • customerservice@vqorthocare.com


The publishers acknowledge that the illustrations used in this publication on pages 2,3,4,5,6,7, 9,10,11,12,13,14,16,17,18,19 are the copyrighted property of VHI Inc.

This operating manual is proprietary and confidential. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording or information storage and retrieval system, without permission in writing from VQ OrthoCare.

Quick Reference Exercise Listing


Shoulder

Shoulder Pendulum (Shoulder Warm-Up)	2
Forward Flexion (Shoulder Flexion)	3
Shoulder Abduction	4
External Rotation	5
Hand Behind Back	6
Side to Side with Wand (Horizontal Abduction/Adduction)	7
Over-the-Door Pulley Range of Motion	8
External Rotation	9
Internal Rotation	10
Shoulder Extension	11
Shoulder Flexion	12
Shoulder Adduction	13
Horizontal Abduction	14
Shoulder Forward Press	15
Diagonal Pull (mid to high)	16
External Diagonal Pull (high to low)	17
External Diagonal Pull (mid to high)	18
Internal Diagonal Pull (mid to high)	19


SHOULDER PENDULUM EXERCISE (SHOULDER WARM UP)

Stand with your knees bent slightly. Bend forward at the waist so that your involved arm is hanging straight down toward the floor. Let your involved arm hang and swing freely in front of you. Shift your weight from side to side allowing your arm to swing slowly from side to side. Relax. Again let your involved arm hang in front of you. Shift your weight so that your arm makes a circular motion (start going in a clockwise direction). Repeat this motion in the opposite (counter clockwise) direction. Do this exercise 3 times daily.


FORWARD FLEXION (SHOULDER FLEXION)

Stand with your knees bent slightly. Grasp the shoulder wand with both hands about shoulder width apart. You can grasp the wand with palms up or down. Start with the wand resting across your hips. While keeping your arms straight bring the wand directly over head until you feel a stretch and hold for 10 seconds. Return slowly to starting position and repeat 10 times. Repeat exercise with palms facing the opposite direction. Do this exercise 3 times daily.


SHOULDER ABDUCTION

Stand with your knees bent slightly. Grasp the wand with both hands. Grasp the wand with your palm up on the injured side and with your palm down on the uninjured side. Start with the wand resting across your hips. Push the wand directly out to the side with the injured side (palm down) until you feel a stretch and hold for 10 seconds. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


EXTERNAL ROTATION

Lie on your back. Grasp the wand with both hands. Grasp the wand with your palm up on the injured side and with your palm down on the uninjured side. Bend both elbows 90 degrees. Keep the elbow of the injured shoulder arm at your side locked at your side. With the arm of the uninjured shoulder, push the wand upward and toward the injured side out as far as you can. Tolerate until you feel a stretch and hold for 10 seconds. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


HAND BEHIND BACK

Stand with your knees bent slightly. Hold towel or wand the way you would hold a towel to dry off your back. The hand on the side of the uninjured shoulder should be at the top of the wand with the hand on the side of the injured shoulder at the bottom of the wand. Use the uninjured arm to slowly pull the injured arm up your back as far as you can tolerate and hold for 10 seconds. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


SIDE TO SIDE WITH WAND **(HORIZONTAL ABDUCTION/ADDUCTION)**

Lie on your back. With palms facing downward, grasp the wand with both hands shoulder width apart. Bring the wand up so your arms are straight out in front of your chest. Leading with the uninjured shoulder, push the wand across your body so that your injured shoulder is moving towards the floor. Then pull back across your body. Do not allow your trunk to twist. You can alternate the grip on the wand with palms facing up then palms facing down. Repeat 10 times. Do this exercise 3 times daily.


OVER-THE-DOOR PULLEY RANGE OF MOTION EXERCISE

Place the over-the-door pulley frame at the top of the doorjamb and close the door to secure. Make sure the pulley is placed on the side of the door that pulls toward you. Grab both handles. Using the uninjured arm, pull downward while you keep the injured arm relaxed. This will raise the affected arm into an extended position. Relax the extended arm and hold the extension for 5-10 seconds and return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily. This exercise can be performed in multiple positions of the shoulder. (For example: You can have your arms directly in front of you, at a 45 degree angle or out to your side at a 90 degree angle).


EXTERNAL ROTATION

Place the exercise band handle between the door and the doorframe about waist level and close the door to secure the handle. Stand sideways with your injured shoulder away from the door. Roll up a towel and place it between your side and upper arm just above your elbow. Grab the handle on the exercise band and make sure your elbow remains bent at 90 degrees as shown in the picture. Rotate your arm outward away from your body keeping your elbow close to your side and return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


INTERNAL ROTATION

Place the exercise band handle between the door and the doorframe about waist level and close the door to secure the handle. Stand sideways with your injured shoulder closest to the door. Place a small towel roll between your side and upper arm just above your elbow. Grab the handle of the exercise band and make sure your elbow remains bent at 90 degrees as shown in the picture. Rotate your arm inward across your body keeping your elbow close to your side and return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


SHOULDER EXTENSION

Place the exercise band handle between the door and the doorframe about waist level and close the door to secure the handle. Stand facing the door with the injured arm extended and grab the exercise band with your palm facing the floor. Keep your arm straight. Pull the handle down and back as far as you can tolerate or just beyond your body as shown and return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


SHOULDER FLEXION

Place the exercise band handle between the door and the doorframe about waist level and close the door to secure the handle. Stand facing away from the door and grab the handle of the exercise band with the injured side. Pull the tubing forward and upward as tolerated or to 90 degrees then return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


SHOULDER ADDUCTION

Place the exercise band handle between the door and the doorframe about waist level and close the door to secure the handle. Stand with your injured shoulder closest to the door and grab the handle of the exercise band with your arm 45 degrees from your side. Pull the band down and towards your side to the buttock. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


HORIZONTAL ABDUCTION

Place the exercise band handle between the door and the doorframe just above waist level and close the door to secure the handle. Stand with the uninjured side facing the door. Grab the handle of the exercise band with the arm on your injured side straight out in front of you. Keep your arm parallel to the floor. Pull your arm across your body and out to the side. Return slowly to the starting position. Do not rotate your body. Repeat 10 times. Do this exercise 3 times daily.


SHOULDER FORWARD PRESS

Place the exercise band handle between the door and the doorframe about shoulder level and close the door to secure the handle. Stand facing away from the door and grasp the handle with your arm at your side and elbow bent as shown. Push the handle out in front of you until your arm is straight. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


DIAGONAL PULL (MID TO HIGH)

Place the exercise band handle between the door and the doorframe about 6" from the top of the door and close the door to secure the handle. Stand sideways with the injured shoulder closest to the door. Grab the handle and step away from the door to take up the slack in the exercise band so that your arm is straight out to your side with your elbow at shoulder height and your palm is facing forward with your thumb up. Pull the handle downward and across your body toward the opposite hip. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.


EXTERNAL DIAGONAL PULL (HIGH TO LOW)

Place the exercise band handle between the door and the doorframe about 6" from the top of the door and close the door to secure the handle. Stand sideways with the injured shoulder farthest away from the door. Grab the handle and step away from the door to take up the slack in the exercise band so that your arm is straight across your body above your shoulder with the palm facing backward. Pull the handle downward and away from your body toward the opposite hip. Return slowly to starting the position. Repeat 10 times. Do this exercise 3 times daily.

EXTERNAL DIAGONAL PULL (MID TO HIGH)

Place the exercise band handle between the door and the doorframe at waist level and close the door to secure the handle. Stand sideways with the injured shoulder farthest away from the door. Grab the handle and step away from the door to take up the slack in the exercise band so that your arm is straight across your body with your palm facing backward. Pull the handle up and across your body and overhead so your palm is now facing forward. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.

INTERNAL DIAGONAL PULL (MID TO HIGH)

Place the exercise band handle between the door and the doorframe at waist level and close the door to secure the handle. Stand sideways with the affected shoulder closest to the door. Grab the handle and step away from the door to take up the slack in the exercise band so that your arm is straight out away from your body with your palm facing forward. Pull the handle up across your body and overhead so your palm faces behind you. Return slowly to the starting position. Repeat 10 times. Do this exercise 3 times daily.

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