
# Home Exercise Kit

Lumbar

User Manual

# VQ OrthoCare Exercise Kit

## Lumbar Kit


### Items included in your kit:

- **Exercise Ball**

(see inflation instructions on back page)

- **Backnobber® II**

- **Air Pump**

- **Exercise Band Kit**

- 2 Handles
- Blue, Green, Red Rubber Bands
- Strap with Loop
- Pouch

- **Towel**

**VQ OrthoCare** is not liable for any misuse or misunderstanding of the exercise product or operating manual. Please call your local representative if any additional assistance is required regarding this product and its operating instructions.

**Caution:** Consult with your therapist or physician before starting any rehabilitation or exercise program. The products in this exercise kit are not to be used by children without adult supervision.

**Caution:** This product contains natural rubber latex which may cause allergic reactions. As with all natural rubber products, keep this product away from grease, oil or excessive heat.

**Disclaimer:** Equipment and supplies illustrated and specified are subject to modification consistent with product availability or other factors.

**VQ OrthoCare**

18011 Mitchell South, Irvine, CA 92614-6007

800.266.6969 • Fax 800.821.8012

[www.vqorthocare.com](http://www.vqorthocare.com) • [customerservice@vqorthocare.com](mailto:customerservice@vqorthocare.com)

The publishers acknowledge that the illustrations used in this publication on pages 2,3,4,5,6,7,10,11,12,13,15,16,17,18 are the copyrighted property of VHI Inc.

Backnobber is a registered trademark of The Pressure Positive Company.

This operating manual is proprietary and confidential. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording or information storage and retrieval system, without permission in writing from VQ OrthoCare.

# Quick Reference Exercise Listing

## *Lumbar*

Pelvic Tilt .....2

Low Back Stretch .....3

Spine Stretch .....4

Upper Back Extension .....5

Prone On Elbows .....6

Exercise Ball Bridge Lying On Floor .....7

Exercise Ball Bridge Single Leg Raise .....8

Abdominal Crunch .....9

Exercise Ball – Reclined Sit-Up .....10

Exercise Ball – Alternate Arm Raise .....11

Exercise Ball – Alternate Leg Raise .....12

Exercise Ball – Alternate Arm and Leg Raise .....13


Exercise Ball – Trunk Rotations .....14

Quadrapped Leg Extension .....15

Quadrapped Alternate Arm and Leg .....16


Exercise Ball – Trunk Extension .....17

Exercise Ball – Prone Hip Extension .....18

## PELVIC TILT

Lie on your back with your knees bent and feet flat on the floor. Flatten your back toward the floor as you tighten your stomach muscles. Relax. Now tilt your pelvis and hips forward creating an arch in your lower back. Breathe out for 5 seconds. Breathe in and relax. Avoid tightening your neck, shoulders and buttocks. Do 15 pelvic tilts. Do this exercise 3 times daily.

## LOWER BACK STRETCH

Lie on your back with your knees bent and feet flat on the floor. Put your hands behind one knee. Pull this knee up toward your chest until a comfortable stretch is felt in the lower back. Hold for 5 seconds. Release this pressure and breathe out slowly. Return the leg to the starting position. Repeat stretch with the alternate leg. Do this stretch 10 times on each leg. Do this exercise 3 times daily.

## SPINE STRETCH


Begin on your hands and knees. Slowly sit back so that you are sitting on your heels. Gently slide your hands forward so that you are stretching your back and supporting your body weight on your outstretched hands. Your forehead and forearms can rest on the floor comfortably. Hold your position and breathe gently in and out for 30 to 60 seconds while you stretch your back. Do this exercise 3 times daily.

## UPPER BACK EXTENSION


Lie face down on the floor. Your forehead may rest on the floor or you can place the folded towel under your forehead. Begin with your arms extended on either side of your body or crossed behind your back. Slowly raise your head and shoulders away from the floor and look up. Come up only as far as your arms comfortably go. Hold for 3 seconds. Slowly lower yourself to a resting position. Relax. Repeat five times. Do this exercise 3 times daily.


## PRONE ON ELBOWS

Lie face down on the floor. Bring your arms up below you to support your weight with your elbows bent at 90 degrees. Push up with your elbows bent and forearms on the floor. Arch your back and raise your torso as high as you can comfortably go. Keep your pelvis on the floor. Hold for 3 seconds. Slowly lower. Relax. Repeat five times. Do this exercise 3 times daily.



## EXERCISE BALL BRIDGE – LYING ON FLOOR

Lie on your back and place the ball under your calves. Tighten your abdominal muscles. Raise your hips off the floor while maintaining abdominal contraction. Do not hold your breath. Do not allow the ball to roll or allow your leg to change position on the ball. Hold for 3-5 seconds. Slowly lower your hips to the floor. Repeat 10-15 times. Do this exercise 3 times daily.


## EXERCISE BALL BRIDGE SINGLE LEG RAISE

Lie on your back and place the ball under your calves. Tighten your abdominal muscles. Raise your hips off the floor while maintaining abdominal contraction. Do not hold your breath. While keeping one leg on the ball, lift one leg about 12 inches towards the ceiling. Hold for 3-5 seconds. Return your leg to ball. Slowly lower your hips to the floor. Repeat 10-15 times for each leg. Alternate legs. Do this exercise 3 times daily.


## ABDOMINAL CRUNCH

Lie on your back and place the ball under your calves. Cross your arms over your chest or place your hands behind your head to make the exercise more difficult. Tighten your abdominal muscles. Do not hold your breath. Do an abdominal crunch so that your shoulder blades lift off the floor. Do not lead with your head and neck. Do not do a full sit-up. Slowly return to the starting position. Repeat 10-15 times. If you place your hands behind your neck do not pull your neck forward. Rest your head in your hands and look towards the ceiling. Do this exercise 3 times daily.


## EXERCISE BALL – RECLINED SIT-UP

Sit on the ball with your back straight and your eyes facing forward with ears lined up over your shoulders. Position your shoulders back as if you are slightly pinching them together and keep your feet apart as wide as your hips. Slowly slide or walk out into a lying position so the ball is under your low back. Keep your knees bent. Cross your arms over your chest or place your hands behind your head to make the exercise more difficult. Tighten your abdominal muscles. Do not hold your breath. Do an abdominal crunch so that your shoulder blades lift off the ball. Do not do a full sit-up. Slowly return to a starting position. Repeat 10-15 times. If you place your hands behind your neck do not pull your neck forward. Rest your head in your hands and look towards the ceiling. Do this exercise 3 times daily.


## EXERCISE BALL— ALTERNATE ARM RAISE

Sit on ball with your head up, eyes facing forward, ears lined up over your shoulders. Keep your shoulders slightly back as if you are slightly pinching them together and keep your feet apart as wide as your hips. Maintain this posture while sitting on the ball by keeping your chest up and maintain the arch in your low back. Tighten your abdominal muscles. Slowly raise one arm overhead and lower to side while maintaining good posture, balance and abdominal contraction. Repeat 10-15 times. Repeat with other arm. Do this exercise 3 times daily.

## **EXERCISE BALL— ALTERNATE LEG RAISE**


Sit on the ball with your head up and eyes facing forward. Have your ears lined up over your shoulders. Keep your shoulders back as if you are slightly pinching them together and keep your feet apart as wide as your hips. Keep proper posture while sitting on the ball by keeping your chest up and maintaining the arch in your low back. Tighten your abdominal muscles. Slowly raise one foot off the floor then lower to floor while maintaining good posture and balance. Alternate with opposite leg. Repeat 10-15 times. Do this exercise 3 times daily.

## EXERCISE BALL— ALTERNATE ARM AND LEG RAISE


Sit on the ball with your head up and eyes facing forward. Have your ears lined up over your shoulders. Keep your shoulders back as if you are slightly pinching them together and keep your feet apart as wide as your hips. Keep proper posture while sitting on the ball by keeping your chest up and maintaining the arch in your low back. Tighten your abdominal muscles. Slowly raise one arm overhead and lift the opposite foot off the floor while maintaining good posture and balance. Alternate leg and opposite arm and repeat 10-15 times on each side. Do this exercise 3 times daily.


## EXERCISE BALL – TRUNK ROTATIONS

Place the exercise band handle between the door and the doorframe about 3 feet from the floor and close the door to secure the band. Sit on the ball with your head up and eyes facing forward. Have your ears lined up over your shoulders. Keep your shoulders back as if you are slightly pinching them together and keep your feet apart as wide as your hips. Using the hand closest to the door grasp the handle of the exercise band and your cross arms over your chest. Maintain proper posture while sitting on the ball by keeping your chest up and maintain the arch in your low back. Tighten your abdominal muscles. Slowly, in a controlled motion, turn your trunk away from the door so that your chest is facing away from the door in a 10 o'clock position for left rotation or a 2 o'clock position for right rotation. Return slowly. Alternate sides. Repeat 10-15 times in each direction. Do this exercise 3 times daily.


## QUADRAPED LEG EXTENSION

Lie face down with the ball under your abdomen with your head and eyes facing the floor while balancing on your toes. Tighten your abdominal muscles. Extend one leg back and straighten knee while maintaining the abdominal contraction. Slowly return to starting position. Alternate each leg. Repeat 10-15 times each leg. Do this exercise 3 times daily.


## QUADRAPED ALTERNATE ARM AND LEG

Lie face down with the ball under your abdomen with your head and eyes facing the floor while balancing on your toes. Tighten your abdominal muscles. Extend one leg back and lift the opposite arm overhead while maintaining the abdominal contraction. Slowly return to the starting position. Alternate each arm and opposite leg. Repeat 10-15 times each leg. Do this exercise 3 times daily.

## EXERCISE BALL–TRUNK EXTENSION

Lie face down with the ball under your abdomen with your head and eyes facing the floor while balancing on your toes. Clasp your hands behind your head. Tighten your abdominal muscles. Raise and extend your trunk so that your lower back contracts and arches slightly. Remember to maintain tight abdominal muscles so you do not over-arch your back. Repeat 10–15 times. Do this exercise 3 times daily.


## EXERCISE BALL – PRONE HIP EXTENSION


Lie face down with the ball under your abdomen and walk your hands forward so that the ball is under your thighs. Tighten your abdominal muscles. Maintain abdominal contraction at all times while the ball is under your thighs. Raise one leg away from the ball toward the ceiling about 6 inches. Do not allow the ball to roll. Return slowly and repeat with the other leg. Do as many as you can or up to 10-15 repetitions for each leg. Do this exercise 3 times daily.

Please refer to the  
Backnobber® II Exercise book  
for specific exercises and instruction.


## Instructions for inflating exercise ball

**Step 1:** Attach short white fill tube to end of air pump.

**Step 2:** Insert tip of white fill tube directly into the exercise ball.

**Step 3:** Using the air pump, inflate the exercise ball until the desired level of firmness is attained. Note: The exercise balls may be "under" inflated to give a soft feel or "over" inflated for a firm feel.

**Step 4:** When inflation is complete, quickly remove the air pump from the exercise ball and insert the long white plug, pushing until flush with the ball surface.


Figure 4a: Long white plug

Figure 4b: Completed inflated exercise ball with plug flush with ball surface.

### Optional Inflation Method

If desired, an alternative air pump such as a bicycle pump or compressor pump can be used to inflate the exercise ball.

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