



OActive™

Off-the-shelf Osteoarthritis Knee Brace

Light weight, low profile, single upright osteoarthritis knee brace with unique angle adjustment to optimize pain relief.

Laterally applied frame and hinge corrects for either medial or lateral compartment and avoids contralateral limb contact during gait. Perfect for bilateral brace users as well as those with limited room between thighs or knees during gait.

- Universal size for the ultimate in easy sizing
- Maximized 3 points of leverage for medial compartment
- Quick release buckles for easy removal
- Lateral application avoids tibial crest pressure issues
- Sili-grip liners and strap pads for optimal suspension
- Available with the BioniCare® Knee System for an unparalleled, synergistic approach to knee osteoarthritis

PRODUCT PART NUMBERS

OACTIVE KNEE BRACE

Patient Height		Left	Right
STANDARD	≤ 6'	08-12006-001	08-12008-001
TALL	≥ 5' 10"	08-12026-001	08-12028-001

*Additional anatomical factors should be considered when selecting appropriate brace length

ADDITIONAL ACCESSORIES

OACTIVE STRAP KIT		06-12000-201		
OACTIVE SUSPENSION WRAP		06-12000-601		
LINER KITS		LEFT 06-12003-101	RIGHT 06-12004-101	
Undersleeve	Size	Part Number		
NEOPRENE	S	VQ-UNDSL3		
	M	VQ-UNDSL4		
	L	VQ-UNDSL5		
	XL	VQ-UNDSL6		
	XXL	VQ-UNDSL7		
	XXXL	VQ-UNDSL8		
	XXXXL	VQ-UNDSL9		
	COTTON/LYCRA	XS	0158000101	
		S	0158000102	
M		0158000103		
L		0158000104		
XL		0158000105		
XXL		0158000106		



OACTIVE KNEE BRACE WITH THE BIONICARE® SYSTEM

TO ORDER CALL 800.652.1135

INDICATIONS

Mild to severe unicompartmental osteoarthritis



CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

