



OActive™

Custom Osteoarthritis Knee Brace

Light weight, low profile, single upright osteoarthritis knee brace with unique angle adjustment to optimize pain relief. Laterally applied frame and hinge corrects for either medial or lateral compartment and avoids contralateral limb contact during gait. Perfect for bilateral brace users as well as those with limited room between thighs or knees during gait. State-of-the-art eCast™ custom measurement system ensures intimate fit for superior comfort and compliance

- Maximized 3 points of leverage for medial compartment
- Quick release buckles for easy removal
- Lateral application avoids tibial crest pressure issues
- Sili-grip liners and strap pads for optimal suspension
- Available with the BioniCare® Knee System for an unparalleled, synergistic approach to knee osteoarthritis

ORDERING INFORMATION

OACTIVE CUSTOM BRACE

Left	Right
08-12996-001	08-12998-001

ADDITIONAL ACCESSORIES

OACTIVE STRAP KIT		06-12000-201		
OACTIVE SUSPENSION WRAP		06-12000-601		
LINER KITS		LEFT 06-12003-101	RIGHT 06-12004-101	
Undersleeve	Size	Part Number		
NEOPRENE	S	VQ-UNDSLV3		
	M	VQ-UNDSLV4		
	L	VQ-UNDSLV5		
	XL	VQ-UNDSLV6		
	XXL	VQ-UNDSLV7		
	XXXL	VQ-UNDSLV8		
	XXXXL	VQ-UNDSLV9		
	COTTON/LYCRA	XS	01580000101	
		S	01580000102	
M		01580000103		
L		01580000104		
XL		01580000105		
XXL		01580000106		



OACTIVE KNEE BRACE WITH THE BIONICARE® SYSTEM

TO ORDER CALL 800.652.1135

Custom Brace Color Chart



Copper



Cherry



Merlot



Forest Green



Lagoon



Deep Lilac



Midnight Blue



Steel Blue



Black



Charcoal



Silver



Champagne



*Chameleon Blue

* There is an upcharge for selecting chameleon blue paint.
Due to normal variations in paint, actual color may vary.

INDICATIONS

Mild to severe unicompartmental osteoarthritis

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

