

Falcon™ Knee Brace

Sleeve

- Comfortprene™ laminate construction inner layer for added comfort
- Polycentric ROM hinges with flexion and extension stops
- Bound edges and Mauser seams for structural integrity
- Popliteal window with spandex material to reduce posterior knee compression
- Breathable spacer fabric around the knee for increased comfort
- Moveable buttress for customized patella control (sleeve only)
- Elasticized proximal strap designed to reduce or eliminate muscle constriction and the resulting distal migration



Wrap



Falcon™ Knee Brace

Product Profile

ROM Hinge



Indications

Mild to Moderate ACL, PCL, MCL, LCL or Combined Ligament Instability, Patellar Subluxation

Warnings

- If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare™, 1.800.452.7993 or consult a physician.
- Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene byproducts.

Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

Sleeve



Wrap



Sizes

	ROM Sleeve	ROM Wrap
Small	07-22015-00-02V	07-22016-00-02V
Medium	07-22015-00-03V	07-22016-00-03V
Large	07-22015-00-04V	07-22016-00-04V
X Large	07-22015-00-05V	07-22016-00-05V
XX Large	07-22015-00-06V	07-22016-00-06V
XXX Large	07-22015-00-07V	07-22016-00-07V



18011 Mitchell South, Irvine, CA 92614
 800.266.6969 Fax 800.821.8012
 www.vqorthocare.com
 customerservice@vqorthocare.com

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.

Comfortprene and Falcon are trademarks and VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005