

Prefabricated Ligament Knee Brace

Whether you brace post-surgically, conservatively or for contact sports, your patients need a knee brace that protects their ligaments. But they want a brace that is comfortable, doesn't slide down their leg, and doesn't impede their motion – even during full activity.



The patented Active Thigh Cuff (ATC) is less restrictive than rigid thigh braces, allowing full flexion of the leg.

A Breakthrough in Suspension

The new prefabricated Hawk™ ligament brace combines the benefits of rigid bracing protection with a new level of suspension and comfort.

Featuring the patented Active Thigh Cuff (ATC) system, the Hawk absorbs soft tissue profile changes during activity. This unique feature allows the brace to exert a consistent skeletal force at the thigh, providing unparalleled suspension and comfort. The strong and lightweight carbon composite shells and aerospace-grade steel hinge provide the strength and protection needed by patients with damaged, deficient or repaired ligaments.

Active Thigh Cuff (ATC) System



These photos demonstrate the action of the ATC system under various conditions (the active straps are in red)



Patient Satisfaction Guarantee

VQ OrthoCare™ guarantees that the Hawk™ brace will provide improved comfort and suspension for your patients. If your patient is not satisfied, please return the brace to VQ OrthoCare within the first 30 days of patient delivery for a full refund.

To obtain prescription forms or order a Hawk brace, please contact your dedicated VQ OrthoCare Representative or call 800.266.6969.



18011 Mitchell South, Irvine, CA 92614
800.266.6969 Fax 800.821.8012
www.vqorthocare.com
customerservice@vqorthocare.com



The new Hawk™ ligament brace features the patented Active Thigh Cuff (ATC) system to dynamically suspend the brace from full flexion to full extension.

Features

- Patented Active Thigh Cuff (ATC)
- Strong and lightweight carbon fiber/plastic composite shells
- Aerospace-grade steel hinge
- Choice of extension stops in 5, 10, 15, 20 or 25 degrees, or full ROM
- Soft suede liners for comfort and suspension
- Selection of right or left, and standard or short length
- Available in Ultra Black finish

Ordering Information

| Hawk – Standard Length | Left | Right |
|----------------------------|-------------|-------------|
| Small | 52405-00-12 | 52405-00-22 |
| Medium | 52405-00-13 | 52405-00-23 |
| Large | 52405-00-14 | 52405-00-24 |
| X-Large | 52405-00-15 | 52405-00-25 |
| XXL | 52405-00-16 | 52405-00-26 |
| Hawk – Short Length | Left | Right |
| Small | 52410-00-12 | 52410-00-22 |
| Medium | 52410-00-13 | 52410-00-23 |
| Large | 52410-00-14 | 52410-00-24 |
| X-Large | 52410-00-15 | 52410-00-25 |
| XXL | 52410-00-16 | 52410-00-26 |
| Accessories | | |
| Undersleeve, XS | 50731-00-01 | |
| Undersleeve, Small | 50731-00-02 | |
| Undersleeve, Medium | 50731-00-03 | |
| Undersleeve, Large | 50731-00-04 | |
| Undersleeve, X-Large | 50731-00-05 | |
| Oversleeve | 50728-00-00 | |
| Patella Cup, Small/Medium | 58003-00-03 | |
| Patella Cup, Large/X-Large | 58003-00-04 | |

Sizing

| Size | Thigh Circumference ¹ | Calf Circumference ² |
|----------|----------------------------------|---------------------------------|
| Small | 15 1/2" - 18" | 12" - 13 3/4" |
| Medium | 18" - 22" | 13 3/4" - 15 1/2" |
| Large | 22" - 24 1/2" | 15" - 16 1/2" |
| X Large | 24 1/2" - 27 1/2" | 16 1/2" - 19 1/4" |
| XX Large | 27 1/2" - 30 1/2" | 19 1/4" - 24 1/4" |

¹ Measured 6" above mid-patella
² Measured 6" below mid-patella

Warnings

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare, 800.452.7993 or consult a physician.

Indications

- ACL, MCL, PCL, LCL and combined instabilities
- Pre- and post-surgery
- Conservative treatment of ligament damage or deficiency
- Prophylactic bracing for contact sports

Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the state in which he/she practices to use or order the use of this device.

Patent number 4,686,969.
Hawk is a trademark and VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005

VQ08051502