

# Portable Hand CPM w/Soft Splint 8091™




HCPM 8091™

## Kinetec® Continuous Passive Motion

- Flexion/extension of fingers
- Aids in preventing post-operative joint stiffness and soft tissue contractures
- Useful for a wide variety of applications
- Anatomically correct motion
- Versatile, comfortable and portable patient therapy
- Patient compliance monitors

# Portable Hand CPM w/Soft Splint 8091™

Kinetec® Continuous Passive Motion



## Product Profile

### Features and Benefits

- Lightweight – .9 lb (397 g) . . . . . Increased patient comfort, improved patient compliance
- Bilateral soft splint . . . . . Improved patient compliance
- Dorsal positioning . . . . . Provides anatomical motion
- Anatomically shaped drive bar . . . . . Further enhances anatomical motion by working finger web spaces
- Non-ski finger wraps . . . . . Help hold fingers securely in place throughout motion
- 15° hyperextension to 270° flexion . . . . . Allows full range of motion
- Up to 15 minutes of pause at extension and flexion limits. . . . . Provides a controlled stretch and rest period
- Composite flexion, intrinsic plus and motion available . . . . . Provides comprehensive motion therapy for treating a wide range of indications
- Digital display of angle, ROM settings . . . . . Convenience
- Timer . . . . . Discreetly monitors patient compliance
- Variable speed and force . . . . . Improves control of therapy
- Safe and reliable

### Specifications

Range of Motion	15° hyperextension to 270° flexion	
Speed	Minimum: 1-2 min	0° - 270° - 0°
	Maximum: 5-30 sec	0° - 270° - 0°
Weight	.9 lb (397 g)	

### Indications

- Open reduction and rigid internal fixation of intra-articular diathyseal and metaphyseal fractures of the phalanges
- Capsulotomy, artholysis and tenolysis for post-traumatic stiffness of MCP and PIP joints
- Crush injuries
- Flexor and extensor tendon synovectomies
- Flexor and extensor tendon tenolysis
- Prosthetic replacement of MCP and PIP joints

### Contraindications/Precautions

- Not to be used in the presence of untreated or uncontrolled infection, unstable fractures or spastic paralysis.
- If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare™, 800.452.7993 or consult a physician.

**CAUTION**  
Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

8091 is a trademark of AbilityOne Corporation.  
Kinetec is a registered trademark of Smith & Nephew SA.  
VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005

(A) Left drive bar	5322-0053 pkg/2	
(B) Right drive bar	5322-0054 pkg/2	
(C) Finger slides	5322-0050 pkg/4	
(D) Lumbrical bill performed	Left	Right
	Small	5322-0043 5322-0046
	Medium	5322-0044 5322-0047
	Large	5322-0045 5322-0048 each
Gloves (not shown)	Left	Right
	Small	5322-0033 5322-0035
	Medium	5322-0034 5322-0036
(E) Finger wraps	5322-0049 pkg/4	
(F) Complete MCP blocking system	5322-0060 each	
(G) MCP block replacement pad kit	5322-0061 pkg/5	
Red clips (not shown)	5322-0055 each	
Wall transformer (not shown)	5322-0056 each	
Carrying case – soft side	5322-0058 each	
Bilateral splint	Small	5322-0068 each
	Large	5322-0069 each
Cover for splint	Small	5322-0070 each
	Large	5322-0071 each

**VQ OrthoCare™**  
18011 Mitchell South, Irvine, CA 92614  
800.266.6969 Fax 800.821.8012  
www.vqorthocare.com  
customerservice@vqorthocare.com