

Thermophore®



Moist Heat Pack

The Thermophore® Automatic Moist Heat Pack™ is a convenient pain-relieving product.



Penetrating moist heat dilates the blood vessels, increasing circulation and creating blood flow. Moist heat therapy is a means to temporarily relieve pain derived from arthritis, bursitis, back aches and muscle spasms.

High temperature and moisture are what make the Thermophore Moist Heat Pack effective.

Thermophore®

Automatic Moist Heat Pack™


Thermophore® uses moisture drawn from the air.



Standard
BAT 075



Medium
BAT 076



Muff
BAT 074

Product Profile

Features and Benefits

- **Effective:** delivers high-temperature moist heat to stimulate circulation, relieve pain and stiffness, and relax muscle spasms or tension.
- **Convenient:** there is no water to boil, no sponge. No scalded fingers, no towels and no need to rely on others to help administer treatments.
- **Efficient:** the Thermophore pack stimulates tissue activity with deep, penetrating heat and temporary pain relief results in a matter of minutes.
- **Safe:** the Thermophore pack is manufactured with a momentary switch which has a feather-light touch for those with arthritic hands. The momentary switch turns off if the patient falls asleep using the product.
- **High Quality:** recognized and prescribed by medical and home care professionals as the best electric fomentation unit on the market.
- **Economical:** simply plug into a 110V outlet. Genuine Thermophore Automatic Moist Heat Packs™ are designed to last for many, many years.
- **Full one-year warranty** for the first year after original purchase, we will repair or replace this product free of charge to the owner if it has a manufacturing defect, or if in use as described by our labels and brochures does not perform as stated.

Adverse Reactions

If used improperly, may cause burns, bleeding or fainting.

Precautions

- Pregnancy
- Areas with poor circulation or poor thermal regulation
- Edema
- Cardiac insufficiency
- Metal
- Open wounds
- Over a topical counterirritant (e.g. FlexAll)

Indications

- Vasodilation
- Reduction of muscle spasms
- Increases local metabolic rate
- Increases tissue extensibility
- Increases pain threshold

Contraindications

- Acute Injury or inflammation
- Recent or potential hemorrhage
- Thrombophlebitis
- Impaired sensation or mentation
- Malignancy

Warnings

This product is to be used under the supervision of a medical professional. This device is not intended for unsupervised public use. If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, consult Patient Care at VQ OrthoCareSM, 1-800-452-7993 or consult a physician.

Ordering Information

Part number	Size	Description	Recommended usage
BAT 075	14" x 27"	Standard	Spine, other large areas of body
BAT 076	14" x 14"	Medium	Knee, shoulder, ankle, small areas
BAT 074		Muff	

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the state in which he/she practices to use or order the use of this device.



18011 Mitchell South, Irvine, CA 92614
800.266.6969 Fax 800.821.8012
www.vqorthocare.com
customerservice@vqorthocare.com

Automatic Moist Heat Pack is a trademark and Thermophore is a registered trademark of Battle Creek Equipment Company
© 2006 VisionQuest Industries, Inc.
VQ0120554REV03