

# Stryker® PainPump®

Comfortable, Continuous, Controlled



## Disposable Syringe Infuser

The Stryker® PainPump® provides continuous, controlled infusion of a local anesthetic directly into the intra-operative site for post-operative pain management.

- Control pain
- Reduce physical therapy
- Consistent delivery
- Reduce patient downtime

# Stryker® PainPump®

## Disposable Syringe Infuser

### Product Profile

#### Features and Benefits

- **120 ml Infuser\***  
Double the volume of most disposable devices. Reduce the need for multiple infusers for one course of therapy.
  - **Burst-Proof Design**  
Non-elastomeric design eliminates the possibility of bursting during filling and use.
  - **Safe and Effective Infusion Monitoring**  
Graduation marks every 5 ml allow patient and clinician to simply and accurately monitor infusion status.
  - **Separate Infuser and Infusion Sets**  
Select the duration of infusion that best meets your therapeutic needs. Simplify stocking requirements.
  - **Designed for Patient Convenience**  
Simple to understand and operate, comfortable to wear, and easy to troubleshoot.
- \*270 ml capacity Infuser also available

#### Ordering Information

Item Number	Description
501-120	5 Day Infusion Kit (0.83 cc/hr/270 ml capacity)
500-120-050	2 Day Infusion Kit (2.08 cc/hr/5.0" fenestrated catheter/120 ml capacity)
500-120-025	2 Day Infusion Kit (2.08 cc/hr/2.5" fenestrated catheter/120 ml capacity)
500-300	Carry Case (for 120 ml pump)
500-375	Carry Case (for 270 ml pump)

#### Indications

The system is indicated for the relief of pain in patients following surgery, by the continuous administration of medication via catheter into the wound site. It is convenient for use by ambulatory patients.

#### Warnings

DISPOSABLE – Destroy after single use. Do not refill or sterilize. Do not overfill device. Follow drug manufacturer's instructions for the medication being used.

#### Contraindications

Not intended for intravenous drug, epidural, infusion, delivery of blood products, TPN or patients with a history of allergic reactions to anesthetics.



18011 Mitchell South, Irvine, CA 92614  
Call 800.266.6969 Fax 800.821.8012  
www.vqorthocare.com  
customerservice@vqorthocare.com

#### CAUTION

Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

Stryker and PainPump are registered trademarks of Stryker Corporation.  
© 2006 VisionQuest Industries, Inc.

VQ0130208REV04