

TROM™

Post-Operative Knee Braces

TROM™ offers state-of-the-art range-of-motion control and locking capabilities. Only TROM allows flexion and extension adjustments at the touch of a button and a quick lock from zero to 90 degrees in 10 degree increments. The ease-of-use and versatility of the TROM can support a variety of protocols following major ligament surgery, meniscal repairs, patella realignment, regenerative chondroplasty and stable femoral fractures.



M2 Telescoping
Cool TROM™

The cool version of the TROM provides the same easy-to-use protected range of motion and locking capabilities as the full foam TROM, but is more ideal for warmer climates. The ease-of-use and versatility of the Cool TROM can support a variety of protocols following major ligament surgery, meniscal repairs, patella realignment, regenerative chondroplasty and stable femoral fractures.



Through an exclusive arrangement with dj Orthopedics, VQ OrthoCare™ offers the **VQ BY DJ** brand of orthopedic products. These products add the value of VQ OrthoCare's exclusive sales, service and patient care network to the DonJoy™ line of quality products.

The FastFit and Rehab TROM are designed as a sized, ready-to-fit brace. They save time in the operating room with a simple one-piece application. The easy-to-use TROM hinge provides controlled range of motion and quick locking capabilities to meet a variety of protocols for major ligament surgeries, meniscal repairs and patella realignments. The FastFit TROM is 24 inches in length and the Rehab TROM is 17 inches in length.



M3 FastFit TROM™
Rehab TROM™

Product Profile



TROM Hinge

Features and Benefits

Telescoping Cool TROM™

- Innovative, easy-to-use hinge
- Quick-lock mechanism from zero to 90 degrees
- Brace can extend to ankle to help prevent migration
- The TROM hinge has extension stops at zero, 10, 20 and 30 degrees
- The TROM hinge has flexion stops at 30, 40, 70, 90 degrees and "open"

FastFit TROM™/Rehab TROM™

- Easy-to-use TROM hinge
- Sized, ready-to-fit
- Quick-release buckles
- Rehab TROM cuffs designed for extended rehabilitation

Ordering Information

Part Number

3-0963-9	Telescoping TROM
3-0842-X	FastFit TROM 24" length
3-0295-X	Rehab TROM 17" length

TROM Conversion Kit (FastFit / Rehab TROM only)

3-0745-X

Sizes Available: S-XL Measurements taken 6" above mid-patella

x=2	S	14" - 18" (36 - 46 cm)	x=3	M	18" - 22" (46 - 56 cm)
x=4	L	22" - 27" (56 - 69 cm)	x=5	XL	27" - 32" (69 - 81 cm)

Accessories

3-0400-9	Adjustable pneumatic condyle pad
3-1770-0	Hand pump

Indications

Supports a variety of protocols following major ligament surgery, meniscal repairs, patella realignment, regenerative chondroplasty and stable femoral fractures.

Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

Warnings

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare, 800.452.7993 or consult a physician.

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.