

VQTM Functional Arm Brace

Standard System

- Quick adjust hand piece allows for supinated, neutral and pronated hand positions
- Lightweight bracing alternative to casting
- Range of motion and fixed position at 30 degree increments
- Allows early passive mobilization and rapid progression to active mobilization



Humeral Brace Attached



VQ™ Functional Arm Brace



Pronation



Supination



Neutral

ROM hinge



Fracture shells



Humeral brace



Product Profile

Ordering Information

Forearm circumference is used as a guide for the standard upper limb brace. Humeral circumference is used as a guide for the humeral brace and humeral brace attachment.

Part #	Sizes	Forearm Circ.	Bicep Circ.
OMN 07-30312-02	Small	7" - 11" (18 - 28 cm)	10" - 13" (25 - 33 cm)
OMN 07-30312-03	Medium	10" - 13" (25 - 33 cm)	12" - 15" (30 - 38 cm)
OMN 07-30312-04	Large	12" - 14" (30 - 35 cm)	14" - 17" (35 - 43 cm)
OMN 07-30312-05	XLarge	13" - 15" (33 - 38 cm)	16" - 19" (41 - 48 cm)

Humeral Attachment

Part #	Sizes
OMN 07-30301-02	Small
OMN 07-30301-03	Medium
OMN 07-30301-04	Large

- The humeral brace can be easily attached to the standard brace on site. Purchase both a standard brace and a humeral brace for on-site modification.
- The humeral brace is designed for use with compression stockinette.
- The immobilization strap allows complete shoulder immobilization.

Indications

- Humeral and forearm fractures
- Post-operative shoulder management
- Biceps repair
- Ulnar nerve transposition

CAUTION

Federal law restricts this device to sale by, or on the order of, a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

Precaution

If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare®, 800.452.7993 or consult a physician.



18011 Mitchell South, Irvine, CA 92614
 Call 800.266.6969 Fax 800.821.8012
 www.vqorthocare.com
 customerservice@vqorthocare.com