


# ULTRASLING® II / ULTRASLING II AB

Both versions of the UltraSling® II provide immobilization for rotator cuff repairs, capsular shifts, Bankhart repairs, glenohumeral dislocation/subluxation, and soft tissue repairs/strains.



UltraSling II

## Shoulder Immobilizer with Abduction Pillow

- Breathable, extra-padded fabric for greater comfort, particularly in longer-term use
- Encourages effective healing by allowing the shoulder and arm to remain in a neutral position
- Helps prevent post-operative internal rotation contractures that can occur in a sling.
- Promotes auxiliary air exchange to reduce risk of secondary infections
- Exercise ball stimulates circulation.
- Easy-open front panel encourages forearm exercises
- Allows arm to be positioned in a variety of positions dependent on post-operative requirements (either more posterior or anterior)
- AB version: large, comfortable pillow gives greater abduction where required (45 or 60 degrees)

# ULTRASLING® II / ULTRASLING II AB



UltraSling® II AB

## Product Profile

### Indications

- Rotator cuff repairs
- Capsular shifts
- Bankhart repairs
- Glenohumeral dislocations/subluxation
- Soft tissue repairs/strains

### Warnings

If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare<sup>SM</sup>, 800.452.7993 or consult a physician.

### Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

### Ordering Information

3-0449-X	Fits Right or Left	UltraSling II
3-0450-X	Fits Right or Left	UltraSling II AB
<b>Sizes available: S-L Measure elbow crease to wrist, fits right or left arm</b>		
x=2	Small	up to 13 1/2" (34.5cm)
x=3	Medium	up to 14 1/2" (37.5cm)
x=4	Large	up to 15 1/2" (39.5cm)



18011 Mitchell South, Irvine, CA 92614  
 Call 800.266.6969 Fax 800.821.8012  
 www.vqorthocare.com  
 customerservice@vqorthocare.com

### CAUTION

Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the state in which he/she practices to use or order the use of this device.

UltraSling is a registered trademark of dj Orthopedics, LLC.  
 VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005

V0033052702