

Static-Pro™ Orthoses and Splinting

Static-Pro™ ROM Elbow Orthosis

- Designed to address range of motion limitations, joint stiffness or soft tissue contractures quickly and efficiently
- Wire-Foam™ cuffs are trimmable and custom moldable for sizing and comfort



Static-Pro™ ROM Wrist Orthosis

- Adjustments are easily achieved with the turn of a knob
- Strong, lightweight construction with a low profile for patient comfort
- Soft, absorbent foam liner
- Simulates manual therapy
- Easy to apply and remove
- Incorporates advanced features technology designed to enhance patient benefit and outcomes



Static-Pro™ Pronation/Supination Orthosis



Static-Pro™ Knee Orthosis



Static-Pro™ Orthoses and Splinting

Product Profile

Intended Use

The DeRoyal® Static-Pro™ orthoses are designed for the treatment of joint stiffness, and established soft tissue contractures, which may have resulted from trauma, immobilization, or neurological disorders. The orthosis should only be used under the supervision of a qualified healthcare

professional. Protocols for use vary depending on the diagnosis and should be determined only after consulting with the patient's physician.

Elbow Orthosis

Description	Product #
Right	4425R
Left	4425L

Wrist Orthosis

Description	Product #
Right	4415R
Left	4415L

Pronation/Supination Orthosis

Description	Product #
Right	4403
Left	4403

Knee Orthosis

Description	Product #
Right	4405R
Left	4405L

Indications

Aids in preventing postoperative joint stiffness and soft-tissue contractures, that may have resulted from trauma, immobilization, or neurological disorders.

Soft Goods (Right)	4405KR
Soft Goods (Left)	4405KL

Precautions

Caution should be used if any of these conditions exist: peripheral vascular disease, healing skin wounds, mild spasticity, heterotrophic ossification, and mild osteoporosis.

If you experience pain, swelling, sensation changes or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at 800.452.7993 or consult a physician.

Contraindications

The Static-Pro orthosis should not be used if any of the following conditions exist: severe spasticity, severe osteoporosis, and/or thrombophlebitis. If pain, numbness, swelling, or skin irritation occurs while wearing the orthosis, the device should be removed. The patient should be re-evaluated and appropriate adjustments made.



18011 Mitchell South, Irvine, CA 92614
800.266.6969
www.vqorthocare.com
customerservice@vqorthocare.com

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by, or on the order of, a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.

Static-Pro and Wire-Foam are trademarks and DeRoyal is a registered mark of DeRoyal Industries, Inc.
© 2007 VisionQuest Industries, Inc.
VQ330529REV04