

Wrist and Thumb Supports



Exolite® Wrist Brace



Exolite's Exo-skeleton and Design Features

The lightweight, low-profile design is 20% less bulky (circumference) than other brands and promotes greater patient compliance. The exo-skeleton is contoured to ensure an anatomically correct fit and is designed to provide smooth, controlled resistance to extreme range of motion. It is pre-shaped into a neutral position for comfortable, limited motion/flexibility. The exo-skeleton provides radial and ulnar styloid relief and more room at the base of the thumb, the palmar crease and 5th metacarpal. The web space is 40% thinner than competing products and is contoured/cushioned with a unique "soft touch" material to reduce irritation. The interior is three-dimensionally molded to create a padded surface that enhances comfort. Our unique sizing system allows the straps to be adjusted to provide an accurate fit and full metacarpophalangeal flexion.

Size and Ordering Information

Right	Left	Wrist Circumference	Palm Circumference	Size
517073	517083	5 1/4" - 6 1/2"	6 3/4" - 8"	S
517075	517085	6 1/2" - 7 3/4"	8" - 9 1/4"	M
517077	517087	7 3/4" - 10"	9 1/4" - 10 1/2"	L



ComfortFORM™ Wrist

Durable lightweight fabric is Lycra® lined for breathability and patient comfort. New and improved preformed aluminum stay helps provide anatomically correct fit and proper support. Ideal for sprains, strains and control of carpal tunnel syndrome symptoms.

Right	Left	Wrist Circumference	Size
3-87282	3-87292	up to 5 1/2"	XS
3-87283	3-87293	5 1/2" - 6 1/2"	S
3-87285	3-87295	6 1/2" - 7 1/2"	M
3-87287	3-87297	7 1/2" - 8 1/2"	L
3-87288	3-87298	8 1/2"+	XL



ComfortFORM Wrist w/Abducted Thumb

Durable lightweight fabric is Lycra lined for breathability and patient comfort. Dual aluminum stays support palmer surface of wrist and extensor surface of thumb. Stays can be adjusted for proper angulation. Contact closure straps. For treatment of sprains, strains, scaphoid injuries, carpal tunnel syndrome, Gamekeeper's Thumb, and deQuervain's syndrome symptoms.

Right	Left	Wrist Circumference	Size
3-87302	3-87312	up to 5 1/2"	XS
3-87303	3-87313	5 1/2" - 6 1/2"	S
3-87305	3-87315	6 1/2" - 7 1/2"	M
3-87307	3-87317	7 1/2" - 8 1/2"	L
3-87308	3-87318	8 1/2"+	XL

Wrist and Thumb Supports



Deluxe Wrist/Forearm Support

10" support features combination loop/lock contact closure and vinyl, foam, tricot construction. Circumferential wrist strap for additional wrist support for sprains, strains, and post-operative rehab of the wrist and distal forearm.

Size and Ordering Information

Right	Left	Wrist Circumference	Size
3-87052	3-87062	up to 5 1/2"	XS
3-87053	3-87063	5 1/2" - 6 1/2"	S
3-87055	3-87065	6 1/2" - 7 1/2"	M
3-87057	3-87067	7 1/2" - 8 1/2"	L
3-87058	3-87068	8 1/2"+	XL



Wrist/Thumb Wrap

Neoprene construction provides compressive wrist support with thumb abduction. One size fits most patients.

Right or Left	Size
3-82600	U



Abducted ThumbSPICA™

Volara™ foam laminate construction with elastic contact closure straps. Effective for deQuervain's syndrome and arthritis.

Right	Left	Size
3-87100	3-87101	U

Warnings

- If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare™, 800.452.7993 or consult a physician.
- Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene byproducts

CAUTION

Applicable Federal, State, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.

Lycra is a registered trademark of E.I. du Pont de Nemours, Volara is a registered trademark of Voltek, Inc., Exolite is a registered trademark of Royce Medical Co. ComfortFORM and ThumbSPICA are trademarks of dj Orthopedics, LLC. VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005



18011 Mitchell South, Irvine, CA 92614
800.266.6969 Fax 800.821.8012
www.vqorthocare.com
customerservice@vqorthocare.com