

VQ™ Gait Walkers

Equalizer® Premium Short-Leg Walker



Equalizer® Premium Low-Top Walker



Features and Benefits

- Contoured strut design allows the gait walker frame to conform to patients' anatomy, ensuring a better fit
- Shock absorbing sole reduces impact of heel strike during ambulation
- Specially engineered rocker bottom is low and wide to promote natural and stable gait
- Wide, anatomically shaped uprights increase stability and comfort
- Straps are pre-threaded in the slots, allowing quick and easy application
- Available in five sizes for men and women
- Available in Black

Ordering Information

Size	Description	Part Number
XSmall	Short Leg	ROY Wo200
XSmall	Low Top	ROY Wo300
Small	Short Leg	ROY Wo400
Small	Low Top	ROY Wo500
Medium	Short Leg	ROY Wo600
Medium	Low Top	ROY Wo700
Large	Short Leg	ROY Wo800
Large	Low Top	ROY Wo900
XLarge	Short Leg	ROY W1000
XLarge	Low Top	ROY W1100

Shoe Size	Mens	Womens
XSmall	2 - 4	3.5 - 5.5
Small	4.5 - 7	6 - 8
Medium	7.5 - 10.5	8.5 - 11
Large	10.5 - 12.5	11.5 - 13.5
XLarge	12.5 +	13.5+

Notice

While every effort has been made to obtain the maximum compatibility of function, strength, durability, and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through the use of this device. This device is a rehabilitation device intended for normal walking only. Do not run, jump or attempt any other high level of activity while wearing this device. Take extra care when ascending or descending steps or stairs.

Indications

- Acute ankle sprains
- Soft-tissue injuries of the lower leg
- Stress fractures of the lower leg
- Stable fractures of the foot and ankle
- Soft-tissue injuries of the lower leg
- Stable delayed-union or non-union fractures of the distal tibia and fibula
- Achilles tendon repairs

Warnings

- This product is to be used under the supervision of a medical professional. This product is not intended for unsupervised public use. If you experience pain, swelling, sensation in, or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at 800.452.7993 or consult a physician.
- Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene byproducts.
- Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity. The determination of when to apply the brace is made solely at the discretion of the treating physician.

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.