

Ankle Supports

Night Splint

The Night Splint can help break the cycle of injury, pain and re-injury by limiting contraction of the soft tissues at night and has also been proven to provide rapid, long-term relief of symptoms. The Night Splint is made from plastic and nylon, with a rubber sole, and has a foam liner that is washable.



Super-Lite A.F.O.

The Super-Lite ankle-foot orthosis is low profile, available in both male and female, left and right size options, and can be trimmed to fit the individual patient. The Super-Lite A.F.O. offers a non-operative approach to the challenging pathology of drop foot.



Ankle Supports

Night Splint

Indications

- Plantar fasciitis
- Achilles tendonitis
- Foot and ankle rehabilitation

Contraindications

- Heel ulcers

Warning

The Night Splint has been provided with a sole; however, it is not intended to be worn during ambulation or weight bearing. The sole is intended to protect the splint only.

Precautions

If you experience any prolonged pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCare™, 800.452.7993 or consult a physician.

Specifications

Size	Small	Medium	Large
Men's shoe	5 - 7 1/2	7 1/2 - 10	> 10
Women's shoe	4 - 6	6 - 9	> 9
Product #	NS2000-10	NS2000-12	NS2000-14

Super-Lite A.F.O.

Indications

- Adult acquired flatfoot secondary to posterior tibial tendon dysfunction
- Lateral ankle instability
- DJD of ankle or rear foot
- Peroneal tendinopathy
- Drop foot

Contraindications

- Severe Charcot arthropathy
- Spastic drop foot
- Spastic and/or uncompensated equinus
- Flexion instability of the knee

Specifications

	Right	Left	Size
Male	3-81150	3-81151	U
Female	3-81160	3-81161	U



18011 Mitchell South, Irvine, CA 92614
Call 800.266.6969 Fax 800.821.8012
www.vqorthocare.com
customerservice@vqorthocare.com

CAUTION

Federal law restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005