

Pneumatic Walker



Pneumatic Walker

- Patented pneumatic system increases stability while decreasing pain and swelling
- Air is easily adjusted to achieve desired compression and custom fit
- Accommodates the swelling patterns that occur throughout the rehabilitation process

Pneumatic Walker



Ordering Information

XSmall	#AW0200
Small	#AW0400
Medium	#AW0600
Large	#AW0800
XLarge	#AW1000

Sizing Information

Size	XSmall	Small	Medium	Large	XLarge
Men's Shoe	2 - 4	4.5 - 7	7.5 - 10.5	10.5 - 12.5	12.5+
Women's Shoe	3.5 - 5.5	6 - 8	8.5 - 11.5	11.5 - 13.5	13.5+

Indications

- Acute ankle sprains
- Soft tissue injuries of the lower leg
- Stress fractures of the lower leg
- Stable fractures of the foot and ankle
- Soft tissue injuries of the lower leg
- Stable delayed union or non-union fractures of the distal tibia and fibula
- Achilles tendon repairs

Contraindications

Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity. The determination of when to apply the brace rests strictly at the discretion of the treating physician.

Warnings

Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene byproducts. If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, contact Patient Care at VQ OrthoCareSM, 1.800.452.7993 or consult a physician.



18011 Mitchell South, Irvine, CA 92614
 800.266.6969 Fax 800.821.8012
 www.vqorthocare.com
 customerservice@vqorthocare.com

CAUTION

Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the state in which he/she practices to use or order the use of this device.

VQ OrthoCare is a service mark of VisionQuest Industries, Inc. © 2005

VQO3405602